


# ViroSurf<sup>zero</sup> ECO-CARTRIDGE

## Surface Disinfectant Alcohol-free


Highly effective surface disinfectant spray designed for use on non-invasive medical device surfaces, including those that are sensitive to alcohol. Contains detergent for added cleaning power. Effective against Bacteria, Yeast, Fungi and Viruses including: HBV, HIV and Rotavirus\*

- The original cartridge based delivery system
- Guaranteed and accurate dosing every time
- Safest system available
- Contact time from only 30 seconds
- Wide effective spectrum
- Wide material compatibility
- Contains detergent







Accurate & risk free dosing



30 second contact time\*



Contains 10 cartridges



makes 5 litres of solution

### KEY FEATURES

Virofex is a concentrated, cartridge based delivery system for alcohol free surface disinfectant.

Designed to remove the need to measure and dispense dangerous chemicals and ultimately save practices both time and money while ensuring that the correct exact dosage is used each and every time - improving safety for both staff and patients.

Each Virofex cartridge contains exactly the right amount of concentrate for one 500ml spray bottle, just add water.

The cartridge is placed inside the neck of the exclusively designed bottle and is automatically activated by screwing the trigger mechanism into place.



Virofex can be used on almost all surfaces within the dental surgery, including non invasive medical devices.

Each cartridge dispenser pack contains 10 x 8ml cartridges which is the equivalent of a 5L/5kg container of ready to use disinfectant.

Packs can easily be wall mounted with the included sticky adhesive pads. The outer surface of the cartridge dispenser pack is coated in Biocote - a unique antimicrobial, silver ion coating, ensuring the dispenser pack can be wall mounted even in designated 'clean' areas.

### COMPOSITION

Aqua, Benzyl-C12-18-Alkyl(dimethyl Chlorides, C12-14-Alkyl[(Ethylphenyl)Methyl]Dimethyl Chlorides, Didecylidimethylammonium Chloride, 2-Propanol (not an active, <1%), Alcohols (not an active) C9-11 Ethoxylated, Dipropylene glycol monomethyl ether, Polyhexamethylene biguanide hydrochloride, Cocamidopropyl Betaine, Tetrasodium EDTA, Dye.

### INSTRUCTIONS FOR USE

Fill the system bottle with a maximum of 500ml clean cold water. Put 1 a cartridge into the neck of the system bottle. Place trigger spray onto system neck by passing the dip tube through the Virofex cartridge and screw down. Gently shake the bottle for a few seconds prior to use. Virofex will remain active for 12 months following activation. Spray surfaces and objects to be disinfected from a distance of approximately 30 cm so that they are completely moistened and let the disinfectant take action over a period of 30 seconds (to achieve full efficacy leave for 5 minutes). Allow surface spray to dry or wipe with clean tissue. For professional use only.

### VARIANTS AND PACKAGING

VIR0001 - Starter pack x 6 per box  
 VIR0002 - Refill pack x 30 per box  
 VIR0003 - Empty 500ml bottle x 10 per box

### EFFECTIVE SPECTRUM AND CONTACT TIMES

Microorganism	Contact Time
<b>Bactericidal</b>	<b>30 sec</b>
Pseudomonas aeruginosa	30 sec
Staphylococcus aureus	30 sec
Enterococcus hirae	30 sec
VRE	1 min
MRSA	1 min
Legionella pneumophila	5 min
<b>Yeasticidal</b>	<b>30 sec</b>
Candida Albicans	30 sec
<b>Viruses</b>	
Influenza B	30 sec
Influenza A H1N1	2 min
Rotavirus A	30 sec
Hepatitis B	1 min
HIV	30 sec
<b>Fungicidal</b>	<b>30 sec</b>
Aspergillus niger	30 sec
<b>Mycobacteria</b>	
Mycobacterium avium	2 min
<b>Spores</b>	
Clostridium difficile (C-diff)	5 min

The following test standards have been achieved: EN 14885, EN 14476, EN 14348, EN 13727, EN 13704, EN 13624.

\*Full effective spectrum available at [www.topdental.co.uk](http://www.topdental.co.uk)

