

Bossklein

www.bossklein.com

Surface Disinfectant Wipe Comparison

All V-WIPE products have been developed for use inside clinical environments on both non invasive medical devices and general high traffic, non-porous surfaces.













| | 200 | Market angle of the control of the c | | W DE 10 | The state of the s | |
|-----------------------------|--|--|---|--|--|--|
| | V-WIPE Classic | V-WIPE Ultra | V-WIPE ZERO Classic | V-WIPE ZERO Ultra | V-WIPE ZERO Universal-XL | V-WIPE ZERO Eco-XL |
| SKU(s) | BOS9783 Tub BOS9784 Refill | BOS1000 Tub BOS1001 Refill | BOS2000 Tub BOS2001 Refill | BOS3000 Tub BOS3001 Refill | BOS4001 | BOS8001 |
| Number of wipes per pack | 200 | 160 | 200 | 160 | 200 | 100 |
| Wipe Size | 200x200mm | 175x155mm | 200x200mm | 175x155mm | 267x210mm | 240x200mm |
| Wipe Thickness/Weight | 18gsm | 40gsm | 18gsm | 40gsm | 34gsm | 40gsm |
| Wipe Material | 100% Polypropylene | 100% Polypropylene | 100% Polypropylene | 100% Polypropylene | 100% Polypropylene | 100% Spunlace Viscose |
| Material Construction | SMS Triple Layer Spunbond/ Meltblown/Spunbond | Meltblown Diamond Pattern | SMS Triple Layer Spunbond/ Meltblown/Spunbond | Meltblown Diamond Pattern | SMS Triple Layer Spunbond/ Meltblown/Spunbond | Hydroentangled parallel lapped natural fibres |
| Total Disinfectant | 620g (3.1g per wipe) | 620g (3.9g per wipe) | 600g (3g per wipe) | 550g (3.4g per wipe) | 1050g (5.25g per wipe) | 543g (5.4g per wipe) |
| Area Coverage (per wipe) | Small to medium | Medium to large | Small to medium | Medium to large | Larger | Larger |
| Active Ingredient/s | Alcohol - Ethanol (63%) | Alcohol - Ethanol (63%) | Quaternary Ammonium Compounds | Quaternary Ammonium Compounds | Quaternary Ammonium Compounds | Quaternary Ammonium Compounds |
| Aroma | Fresh Lemon | Fresh Lemon | Subtle Apple | Subtle Apple | Subtle Apple | Subtle Apple |
| EN14476 Tested | Yes | Yes | Yes | Yes | Yes | Yes |
| EN16615 Tested | Yes | Yes | Yes | Yes | Yes | Yes |
| Contact Time | 60 sec | 60 sec | 60 sec | 60 sec | 60 sec | 2 mins |
| 2 in 1: Cleans & Disinfects | No | No | Yes | Yes | Yes | Yes |
| Summary | Quick drying Streak free Most economical wipe Only for use on alcohol resistant surfaces | Quick drying Streak free Premium feel wipe Only for use on alcohol resistant surfaces | All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Most economical wipe | All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Premium feel wipe | All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Largest wipe Strongest wipe | Biodegradable wipe Renewable resources All non porous surfaces Cleans & disinfects Alcohol Free Low streaking |