



Infection Control

This useful guide of a typical dental surgery demonstrates the ideal areas to use our key products, Virofex has **all** bases covered when it comes to hygiene and disinfection.



ViroTorq

Rotary Instrument Disinfectant

- Contact time from only 5 mins*
- Ready to use solution
- For burs, diamonds and polishers
- For use inside bur bath
- Aldehyde and phenol free
- Wide material compatibility



ViroSonic

Ultrasonic Cleaner Concentrate

- 3% concentrate (30ml/1L)
- For pre-cleaning of instruments
- Use inside Ultrasonic bath
- Removes stubborn debris
- Enzymatic
- Wide material compatibility



Quantum

Instrument Disinfectant Colour Change Powder

- 2% concentrate (20g/1L)
- Highly effective
- Economical
- Enzymatic
- EN14476 virus tested
- Effective against Coronavirus



ViroTech

Instrument Disinfectant Concentrate

- From 2% concentrate (20ml/1L)*
- EN14476 virus tested
- Wide material compatibility
- Active against TB
- Effective against Coronavirus
- Mild and fresh aroma



ViroSurf

Surface Disinfectant Alcohol based

- Contact time of only 60 seconds*
- EN14476 virus tested
- Effective against Coronavirus
- Lemon aroma
- Ready to use



ViroSurf wipes

Surface Disinfectant Wipes Alcohol based

- Contact time of only 60 seconds*
- EN14476 virus tested
- Effective against Coronavirus
- Lemon aroma
- Ultra Soft (160 premium wipes)
- Extra Value (250 economy wipes)



ViroGel

Hand Disinfectant Gel

- 30 second contact time
- EN14476 Virus tested
- Effective against Coronavirus
- Clear gel with natural aroma
- Leaves no sticky residue



ViroRub

Hand Disinfectant Rub

- 30 sec contact time
- EN14476 Virus tested
- Effective against Coronavirus
- Clear spray natural aroma
- Quick drying



ViroVent³⁶⁵

Aspirator Disinfectant Daily Concentrate

- Contact time from only 15 mins
- 2% concentrate (20ml/1L)
- Designed for daily use
- Low-foaming and non-abrasive
- Wide effective spectrum
- Chloride, aldehyde & phenol free



ViroVent⁵²

Aspirator Disinfectant Weekly Concentrate

- Contact time of 30 minutes
- 5% concentrate (50ml/1L)
- Designed for weekly use
- Low foaming and non abrasive
- Fully bactericidal
- Chloride, aldehyde & phenol free



ViroImpress

Impression Disinfectant Powder


- Fast acting - just 60 seconds
- 2% concentrate (20g/1L)
- Powder based disinfectant
- Universal compatibility
- EN14476 tested
- Effective against Coronavirus



ViroLine

Dental Waterline Cleaner Concentrate

- 2% concentrate (20ml/1L)
- Dissolves protein
- Removes biofilm
- Keeps water flowing freely
- Aldehyde-free



ViroSurf^{zero}

ECO-CARTRIDGE Surface Disinfectant Alcohol free

- Cartridge based system
- No measuring required
- EN14476 virus tested
- Effective against Coronavirus
- 2 in 1 - cleans and disinfects
- Safe to use on all surfaces



ViroSurf^{zero}

Surface Disinfectant Alcohol free

- Contact time of only 30 secs*
- EN14476 virus tested
- Effective against Coronavirus
- Mint aroma
- 2 in 1 - Cleans and disinfects
- Safe to use on all surfaces



ViroSurf^{zero}

wipes Surface Disinfectant Wipes Alcohol free

- Contact time of only 30 seconds*
- EN14476 virus tested
- Effective against Coronavirus
- Mint aroma
- 2 in 1 - Cleans and disinfects
- Safe to use on all surfaces
- Extra Value (250 economy wipes)

instruments
HAND & ROTARY

surfaces
ALCOHOL-BASED ALCOHOL-FREE

special areas
ASPIRATORS, WATER LINES & IMPRESSIONS

hand hygiene
GEL & RUB